



Xpac – Powerful independence

Xpac by Breas is an extended-life backup battery with a difference. The modular design contains two high powered lithium ion batteries, which can be easily removed and replaced by the user. The Xpac intuitive LED display shows the health of each battery so you know when it is time to renew the batteries.



LONGER RUN TIMES

- Double the capacity of similar external batteries, containing 2 x 99.4 Wh batteries

LOWER COST OF OWNERSHIP

- Removable batteries which can be replaced when they reach the end of life

INTELLIGENT AND INTUITIVE LED DISPLAY

- The battery management system allows users to check the state of charge/health

QUICK INDEPENDENT CHARGING SYSTEM

- Charges both batteries from empty in <4 hours

MOBILITY

- Option of shoulder carry bag which also protects the battery when placed on/under wheelchair
- Certified for use on aircrafts (remove one battery)

Breas Xpac Technical Specifications

Approximate Run Times				
Device	Vivo 30/40	Vivo 50/55/60/65	Vivo 45/45LS Nippy 4/4+	Vivo 1-2-3
Run Time	10 hours*	12 hours*	18 hours*	36 hours*
Recharge Time	4 hours (both batteries charge simultaneously)			
* IPAP 20 cm/EPAP 8 cm, Rate 15 bpm				

Technical Information	
Battery Type	Lithium Ion
Battery Capacity	99.4 Wh (per battery)
Output	19v -34.8v
Power Supply	Input 100 – 240v Output 18v DC, 60w Max
Weight (Dual Battery)	1.5 kg (3.3 lb)
Dimensions L x H x D	213 x 58 x 94 mm (8.4 x 2.3 x 2.7 inches)
Life Span	Battery 2 – 3 years (with daily use)

Order Codes				
	USA	EU	UK	JP
Xpac Single	007994	007995	007993	007996
Xpac Dual	007998	007999	007997	008000
Carry Bag	007932			

Spare Parts	
Power Supply	007727
Spare Lids	009021
Main Chassis	007806

Cables



Xpac by Breas is a trademark of Breas Medical AB. © 2022 Breas Medical – All rights reserved. Breas Medical reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information. Breas and the Breas logo are trademarks of Breas Medical AB. COR-MAR-MAR-007851-Rev3